


Anleitung Blob (scroll down for english version)

1. **Tipp:** Teile das Sklerotium in 2 Teile, falls mit dem ersten Blob etwas schief geht, hast du noch einmal Ersatz.



2. Befeuchte das Sklerotium mit Wasser (am besten destilliert) und lege es mit dem Blob nach oben in eine Schale.



3. Der Blob mag kein helles Licht, vor allem keine Sonne!
Decke die Schale mit einem dicken Tuch oder Kleidungsstück ab.
4. Nach einigen Stunden sollte er beginnen sich auszubreiten.



5. Lege in Richtung seiner Ausbreitung ein oder mehrere Haferflocken, um den Blob zu füttern.



allgemeine Hinweise:

- Der Blob mag am liebsten $25\text{ °C} \pm 5\text{ °C}$.
- Er muss immer leicht feucht gehalten werden, sonst entwickelt er sich zum Sklerotium.
- Wenn die alten Haferflocken schleimig und verdorben aussehen, sollte ein Teil des frischen Blobs in eine neue Schale gesetzt werden.
- Am besten füllt man in die Schale eine 3 mm hohe **AgarAgar**-Schicht, darauf wächst er besser und man kann ihn einfacher in eine andere Schale umsetzen. *Dazu 0,2 g Agar + 100 ml destilliertes Wasser 3 min lang köcheln lassen und noch heiß in die Schale gießen.*

Viel Spaß mit dem Blob!




Instructions for the blob

1. **Tip:** Divide the sclerotium into 2 parts. In case something is going wrong with the first blob, you have one more sclerotium.



2. Moisten the sclerotium with water (preferably distilled) and put it with the blob facing up in a petridish or bowl.



3. The blob doesn't like bright light, especially not the sun! Cover the dish with a thick cloth or piece of clothing.
4. After a few hours it should begin to spread.



5. Put one or more oatmeal in the direction of its expansion to feed the blob.



General information:

- The blob likes $25\text{ }^{\circ}\text{C} \pm 5\text{ }^{\circ}\text{C}$ best.
- It must always be kept slightly moist, otherwise it will develop into a sclerotium.
- If the old oatmeal looks slimy and spoiled, some of the fresh blob should be placed in a new dish.
- It is best to fill a 3 mm high **layer of agar agar** into the dish, it grows better on it and it's easier to move it to another dish.
To do this, simmer 0.2 g agar + 100 ml distilled water for 3 minutes and pour it into the dish while it is still hot.

Have fun with the Blob!

